RetroNectin® Pro

(Recombinant Human Fibronectin Fragment)

Code No. T101A

Size: 0.5 mg / 0.5 ml

Code No. T101B

Size: 2.5 mg / 2.5 ml

Description:

RetroNectin Pro (Recombinant Human Fibronectin Fragment) is a recombinant human fibronectin fragment (rFN-CH-296) and composed of three functional domains: the cell-binding domain (C-domain), heparin-binding domain (H-domain), and CS-1 sequence (Figure). The fragment enhances retroviral/lentiviral-mediated gene transduction by aiding the co-localization of target cells and virions. Specifically, virus particles bind RetroNectin via interaction with the H-domain, and target cells bind mainly through the interaction of cell surface integrin receptor VLA-5 and VLA-4 with the fibronectin C-domain and CS-1 site, respectively. By facilitating close proximity, RetroNectin can enhance retroviral-mediated gene transfer to target cells. In addition, the reagent enhances the proliferation of T lymphocytes. T cells are conventionally expanded in the presence of Interleukin-2 (IL-2) by stimulation with anti-CD3 antibody. The addition of RetroNectin in this stimulation step dramatically increases the efficiency of T cell expansion. Moreover, T cells expanded with this method contain a high proportion of less-differentiated T cells (naive T cells) that have the ability to differentiate into cytotoxic T cells *in vivo*.

Through the optimization of the manufacturing process, RetroNectin Pro is suitable for research and development in the field of cell and gene therapy such as CAR-T cell therapy. For more information on the RetroNectin method, please refer to technical literature and product information available online.


Figure. The hypothesized mechanism of RetroNectin-mediated transduction. The cell binds to the CS-1 site via VLA-4, and to the C-domain via VLA-5. The viral particle can bind to the H-domain of RetroNectin. These interactions increase the localized concentrations of cells and viral particles, an effect that is thought to enhance gene transduction.

Storage: -20°C

Molecular Mass: 62,617 Da (amino acid sequence)

Source: *E. coli* expressing human fibronectin fragment CH-296

Purity: >90% by HPLC

Form: RetroNectin solution containing 12.5 mM sodium citrate (pH 6.2) and 1.25% sucrose

Packaging: RetroNectin is provided as a sterile 1 mg/ml solution.

Applications: To enhance the efficiency of gene transfer into mammalian cells and enhance the proliferation of T lymphocytes.

Materials Required but not Provided:

[Reagents]

- Sterile PBS (-)
- HBSS/Hepes (Hank's Balanced Salt Solution supplemented with 2.5% (v/v) 1 M Hepes)
- 2% bovine serum albumin (BSA Fraction V)/PBS Solution

[Equipment]

- Non-treated tissue culture plates or dishes
- Electric pipetter
- Pipetter
- Sterile pipettes
- Sterile tips with filters
- Safety cabinet or clean workstation

Precautions for Use:

- Work under sterile conditions.
- Freezing and thawing can be repeated up to 10 times.
- Do not mix the solution vigorously. Do not vortex.
- To avoid loss of RetroNectin fragment, do not filter-sterilize RetroNectin solution diluted with PBS.

Protocols:

[Preparation of RetroNectin-Coated Plates]

Coat a plate using 20 - 100 $\mu\text{g/ml}$ RetroNectin with a volume corresponding to 4 - 20 $\mu\text{g/cm}^2$ plate area.

1. Prior to coating, prepare a RetroNectin solution (20 - 100 $\mu\text{g/ml}^*$) by diluting with sterile PBS.
2. Dispense an appropriate volume^{*2} of sterile RetroNectin solution into each plate and allow the plate to stand for 2 hours at room temperature or at 4°C overnight.
Note: Non-treated, cell culture-grade tissue culture plates or dishes should be used in this step.
3. Remove the RetroNectin solution and then block with an appropriate volume^{*2} of sterile 2% bovine serum albumin (BSA, Fraction V) in PBS. Allow the plate to stand at room temperature for 30 minutes^{*3}.
4. Remove the BSA solution, and wash the plate once with an appropriate volume of HBSS/Hepes or PBS. After removing the wash solution, the plate is ready for use. The RetroNectin-coated plate can be sealed with Parafilm and stored at 4°C for up to one week.

^{*1} Example: calculating the amount of RetroNectin reagent:

When 2.25 ml of RetroNectin solution at a concentration of 20 $\mu\text{g/ml}$ is placed into a 35 mm diameter dish (9 cm^2), the concentration used for coating is 5 $\mu\text{g/cm}^2$.

^{*2} 0.5 ml for each well of a 24-well plate or 2 ml for each well of a 6-well plate.

^{*3} If the plates or the dishes are to be used immediately, the blocking step with 2% BSA/PBS solution can be omitted. In this case, repeat the washing with PBS or HBSS/Hepes in Step 4. twice.

Quality Statement:

This product does not contain any human- or animal-derived materials in the final formulation.

References:

- 1) Kimizuka F, Taguchi Y, Ohdate Y, Kawase Y, Shimojo T, Hashino K, Kato I, Sekiguchi K, and Titani K. Production and characterization of functional domains of human fibronectin expressed in *Escherichia coli*. *J Biochem.* (1991) **110**: 284-291.
- 2) Hanenberg H, Xiao XL, Dilloo D, Hashino K, Kato I, and Williams DA. Colocalization of retrovirus and target cells on specific fibronectin fragments increases genetic transduction of mammalian cells. *Nat Med.* (1996) **2**: 876-882.
- 3) Hanenberg H, Hashino K, Konishi H, Hock RA, Kato I, and Williams DA. Optimization of fibronectin-assisted retroviral gene transfer into human CD34⁺ hematopoietic cells. *Hum Gene Ther.* (1997) **8**: 2193-2206.
- 4) Pollok KE, Hanenberg H, Noblitt TW, Schroeder WL, Kato I, Emanuel D, and Williams DA. High-efficiency gene transfer into normal and adenosine deaminase-deficient T lymphocytes is mediated by transduction on recombinant fibronectin fragments. *J Virol.* (1998) **72**: 4882-4892.
- 5) Chono H, Yoshioka H, Ueno M, and Kato I. Removal of inhibitory substance with recombinant fibronectin-CH-296 plates enhances the retroviral transduction efficiency of CD34⁺CD38⁻ bone marrow cells. *J Biochem.* (2001) **130**: 331-334.
- 6) Hosoi H, Ikeda H, Imai N, Amaie C, Wang L, Orito Y, Yamane M, Ueno H, Ideno M, Nukaya I, Enoki T, Mineno J, Takesako K, Hirano S, and Shiku H. Stimulation through very late antigen-4 and -5 improves the multifunctionality and memory formation of CD8⁺ T cells. *Eur J Immunol.* (2014) **44**: 1747-58.

Related Products:

Anti-CD3 mAb GMP grade (Cat. #T210)

LymphoONE® T-Cell Expansion Xeno-Free Medium, 1L Bottle (Cat. #WK552S)*

* Not available in all geographic locations. Check for availability in your area.

RetroNectin and LymphoONE are registered trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

RetroNectin® Pro

(Recombinant Human Fibronectin Fragment)

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