RetroNectin® GMP grade

Code No. T202

Size: 2.5 mg / 2.5 ml

Description:

RetroNectin, Recombinant Human Fibronectin Fragment CH-296, is a chimeric peptide produced in *E. coli*. This peptide consists of three functional domains: a cell-binding (C-domain), a heparin-binding (H-domain), and a CS-1 site.

RetroNectin GMP grade is manufactured as quality-assured product, according to relevant PIC/S guidelines for Good Manufacturing Practice (GMP).

The RetroNectin method has been widely employed for retroviral and lentiviral vector-mediated gene transfer into mammalian cells. For more details about the RetroNectin method, please refer to technical literature and product information available online for research-grade RetroNectin (Recombinant Human Fibronectin Fragment) (Cat. #T100A/B).

Storage: Store at or below -30°C

Molecular Mass: 62,617 Da (amino acid sequence, 574 amino acids)

Source: *E. coli* expressing human fibronectin fragment CH-296

Purity: >90% by HPLC

Form: 12.5 mM sodium citrate, pH 6.2 and 1.25% sucrose

Packaging: 2.5 mg of protein/2.5 ml in a clear glass 5 ml vial with a slotted stopper and a flip-cap seal

Applications:

To enhance the efficiency of gene transfer into mammalian cells.

Instructions For Use:

[Filtration]

- 1) Materials required
 - Sterile syringe and needle
 - 0.22 μ m Filter (e.g., MERCK MILLIPORE MILLEX-GV, SLGV 033 RS)

- 2) Protocol
Filter the thawed solution (1.0 mg/ml) through a 0.22 μ m Filter.

Note: To avoid protein loss, never filter diluted solution. Only stock solution (1.0 mg protein/ml) should be filtered. Additionally, use low protein binding filters.

[Dilution and RetroNectin Coating]

- 1) Materials required
 - PBS(-) or ACD-A
 - Sterilized pipettes and filter tips

- 2) Protocol
Coat the surface of containers, such as culture dishes, petri dishes, flasks, and bags with RetroNectin Solution, which concentration is 20 - 100 μ g/ml. Also the proper amount of RetroNectin Solution is from 4 - 20 μ g/cm².

1. Before coating, dilute RetroNectin Solution with sterilized phosphate buffered saline (PBS) or Acid Citrate Dextrose Formula A (ACD-A) to adjust its appropriate concentration (ranging 20 - 100 μ g protein/ml).
2. Dispense the appropriate volume of the sterilized RetroNectin solution into containers, then let stand for at least 0.5 hours in room temperature, or overnight at 4°C.
3. Remove the RetroNectin solution, and then add wash solution (e.g. PBS).
4. Remove the wash solution. The surface of the container has been coated, and it is ready to use.

Precautions:

- Work under sterile conditions.
- The thawed product should be used within 12 hours at room temperature.
- Once thawed, use the necessary quantity of thawed solution and discard any unused solution.

Quality Statement:

1. RetroNectin GMP grade is manufactured according to relevant PIC/S GMP guidelines.
2. No animal- or human- derived components are used in the manufacture of this product.
3. Takara Bio has submitted a Drug Master File (DMF application number 18898) to the Food and Drug Administration for this product.

References:

- 1) Kimizuka F, Taguchi Y, Ohdate Y, Kawase Y, Shimojo T, Hashino K, Kato I, Sekiguchi K, and Titani K. Production and characterization of functional domains of human fibronectin expressed in *Escherichia coli*. *J Biochem.* (1991) **110**: 284-291.
- 2) Hanenberg H, Xiao XL, Dilloo D, Hashino K, Kato I, and Williams DA. Colocalization of retrovirus and target cells on specific fibronectin fragments increases genetic transduction of mammalian cells. *Nat Med.* (1996) **2**: 876-882.
- 3) Hanenberg H, Hashino K, Konishi H, Hock RA, Kato I, and Williams DA. Optimization of fibronectin-assisted retroviral gene transfer into human CD34⁺ hematopoietic cells. *Hum Gene Ther.* (1997) **8**: 2193-2206.
- 4) Pollok KE, Hanenberg H, Noblitt TW, Schroeder WL, Kato I, Emanuel D, and Williams DA. High-efficiency gene transfer into normal and adenosine deaminase-deficient T lymphocytes is mediated by transduction on recombinant fibronectin fragments. *J Virol.* (1998) **72**: 4882-4892.
- 5) Chono H, Yoshioka H, Ueno M, and Kato I. Removal of inhibitory substance with recombinant fibronectin-CH-296 plates enhances the retroviral transduction efficiency of CD34⁺CD38⁻ bone marrow cells. *J Biochem.* (2001) **130**: 331-334.

Related Products:

[Research grade RetroNectin, not for clinical trial]

RetroNectin® (Recombinant Human Fibronectin Fragment) 0.5 mg (Cat. #T100A)
RetroNectin® (Recombinant Human Fibronectin Fragment) 2.5 mg (Cat. #T100B)

[Other]

RetroNectin® Dish (RetroNectin Pre-coated Dish, 35 mm ϕ) (Cat. #T110A)
RetroNectin® EIA Kit (Cat. #MK140)
Wash and Stop Solution for ELISA without Sulfuric Acid (Cat. #MK021)
Anti-CD3 mAb GMP grade (Cat. #T210)
LymphoONE™ T-Cell Expansion Xeno-Free Medium, 1L Bottle (Cat. #WK552S)

RetroNectin is a registered trademark of Takara Bio Inc.
LymphoONE is a trademark of Takara Bio Inc.

Note

This product is to be used for research use or non-commercial manufacturing for clinical research. This product is not intended for humans or animals in-vivo applications. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or e-mail: bio-licensing@takara-bio.co.jp. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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