PrimeCap™ T7 RNA Polymerase (low dsRNA) GMP grade

Code No. T2422S001 **Size:** 200 KU
 Conc.: 200 U/μl

Code No. T2422S010 **Size:** 2,000 KU
 Conc.: 200 U/μl

Code No. T2422S100 **Size:** 20,000 KU
 Conc.: 200 U/μl

Description:

This enzyme is a modified form of T7 RNA polymerase that combines the features of cap analog-dependent RNA synthesis, low dsRNA generation, and heat resistance (~52°C). *In vitro* transcription (IVT) using this enzyme, a cap analog, and 10X IVT Reaction Buffer GMP grade (Cat. #T2401S010/S100/S500) can produce a large quantity of high-quality mRNA with high-capping efficiency while significantly reducing the generation of immunogenic dsRNA.

Quality Statement:

1. This product is manufactured in accordance with relevant GMP guidelines. For more information, please refer to our website.
2. No animal- or human-derived components are used in the manufacture of this product.
3. This product uses primary materials that do not contain β-lactam compounds.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Storage: -20 ± 5°C

Source:

Escherichia coli carrying a plasmid containing the gene for phage T7 RNA polymerase variant

Properties:

Molecular mass : approx. 99.8 kDa
Cofactor : Mg²⁺

Unit Definition:

One unit is defined as the amount of enzyme that generates 0.25 μg of 1.9 kb FLuc RNA in 1 hour at 37°C.

Reaction Mixture for Unit Definition:

| | |
|--------------|-------------------------------|
| 1X | IVT Reaction Buffer GMP grade |
| 10 mM | ATP, CTP, GTP, and UTP |
| 0.9 μg/20 μl | Linearized FLuc plasmid DNA |

Applications:

1. Synthesis of capped mRNA using a cap analog (Co-transcriptional capping)
 2. Synthesis of uncapped mRNA for enzymatic capping
- Note: The RNA yield may be reduced by 5 - 30% compared to when using cap analogs.

Precautions for Use:

1. Do not mix the enzyme vigorously.
2. Avoid introducing RNase into the reaction. The presence of RNase can cause RNA fragmentation and lead to significantly reduced RNA yield from the reaction.
3. In order to synthesize a uniform length of RNA, linear dsDNA such as cleaved plasmid and PCR product is used as the template. We recommend that the 3' end of the template DNA should be 5'-protruding end or blunt end to avoid unwanted products.
4. The 10X IVT Reaction Buffer GMP grade contains spermidine, which forms a complex with nucleic acids that may cause insoluble precipitate. To minimize this effect, template DNA should be added just prior to adding the enzymes.

Application Example

(Synthesis of 1.9 kb capped mRNA with CleanCap Reagent AG):

| | |
|--|------------|
| Sterile purified water | x μl |
| 10X IVT Reaction Buffer GMP grade | 2 μl |
| ATP, CTP, GTP, UTP*1 | each 10 mM |
| CleanCap Reagent AG*2 | 4 mM |
| Template DNA | 0.5 - 2 μg |
| RNase inhibitor | 20 U |
| Pyrophosphatase (inorganic) | 0.1 U |
| PrimeCap T7 RNA Polymerase (low dsRNA) GMP grade | 200 U |
| Total | 20 μl |

Incubate at 37°C for 1 - 2 hrs.

*1 When using a modified NTP, replace the corresponding NTP with an equivalent amount of the modified version.

*2 CleanCap Reagent AG (TriLink BioTechnologies: Code. N-7113-1/5/10 etc.)

Related Products:

10X IVT Reaction Buffer GMP grade (Cat. #T2401S010/S100/S500)
T7 RNA Polymerase GMP grade (Cat. #T2411S001/S010/S100)

PrimeCap is a trademark of Takara Bio Inc.

Note

This product is not intended for humans or animals *in vivo* applications. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

PrimeCap™ T7 RNA Polymerase (low dsRNA) GMP grade

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