# Pyrophosphatase (inorganic) GMP grade

Code No. T2441S001      Size:            1,000 U  
   Conc.:            1 U/ $\mu$ l

## Description:

This enzyme exhibits activity comparable to that of Pyrophosphatase (inorganic), HQ (Cat. #2451A) and hydrolyzes inorganic pyrophosphate (PPi) to orthophosphate. By hydrolyzing PPi generated during nucleotide polymerization, it drives enzymatic nucleic acid synthesis forward. When used with T7 RNA Polymerase GMP grade (Cat. #T2411S001/S010/S100) or PrimeCap® T7 RNA Polymerase (low dsRNA) GMP grade (Cat. #T2422S001/S010/S100) in an *in vitro* transcription (IVT) reaction, this enzyme enhances RNA yield.

## Quality Statement:

1. This product is manufactured in accordance with relevant GMP guidelines. For more information, please refer to our website.
2. No animal- or human-derived components are used in the manufacture of this product.
3. This product uses primary materials that do not contain  $\beta$ -lactam compounds.

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Storage:      -20  $\pm$  5°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for *Saccharomyces cerevisiae* (yeast) inorganic pyrophosphatase

## Properties:

Molecular mass : approx. 66 kDa (homodimer)  
Cofactor            : Mg<sup>2+</sup>

## Unit Definition:

One unit is defined as the amount of enzyme that generates 1  $\mu$  mol of phosphate from inorganic pyrophosphate per minute at 25°C.

## Reaction Mixture for Unit Definition:

100 mM	Tris-HCl, pH 7.2
2 mM	MgCl <sub>2</sub>
2 mM	inorganic pyrophosphate (PPi)

## Applications:

1. RNA yield enhancement in *in vitro* transcription reaction
2. Promotion of DNA amplification reaction

## Precautions for Use:

Do not mix the enzyme vigorously.

## References:

- 1) Cooperman B S. *Methods Enzymol.* (1982) **87**: 526-548.
- 2) Cunningham P R and Ofengand J. *Biotechniques.* (1990) **9**: 713-714.
- 3) Dean B F., *et al. Genome Res.* (2001) **11**: 1095-1099.

## Related Products:

PrimeCap® T7 RNA Polymerase (low dsRNA) GMP grade  
(Cat. #T2422S001/S010/S100)  
T7 RNA Polymerase GMP grade (Cat. #T2411S001/S010/S100)  
10X IVT Reaction Buffer GMP grade (Cat. #T2401S010/S100/S500)  
Recombinant RNase Inhibitor GMP grade (Cat. #T2431S001/S010/S100)  
BspQ I GMP grade (Cat. #T2451S001/S010/S100)  
10X BspQ I Buffer GMP grade (Cat. #T2461S010/S100)

PrimeCap is a registered trademark of Takara Bio Inc.

## Note

This product is not intended for humans or animals *in vivo* applications. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# Pyrophosphatase (inorganic) GMP grade

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