



## TECH NOTE

# *In Vitro* Cleavage Efficiency of sgRNAs Correlates with Functional Genome Editing in Target Cells

**A novel *in vitro* assay to test sgRNA cleavage efficiency**

Screen various sgRNAs to determine the most effective sgRNAs prior to delivering to your cells >>

**Accurate prediction of sgRNA cleavage efficiency**

sgRNA cleavage efficiency predicted *in vitro* correlates with *in vivo* cleavage as assessed by both a nuclease assay and functional analysis >>

---


## Overview

In **CRISPR/Cas9 genome editing**, targeting the Cas9 nuclease to a specific genomic locus is solely mediated by a user-defined sgRNA. Currently available web-based tools for sgRNA design will return a variety of candidate sgRNAs for a single gene target. Despite these *in silico* predictions, not every sgRNA will exhibit equivalent cleavage efficiency. Given this inconsistency, it is necessary to screen multiple sgRNAs to identify the most effective one.


## An *In Vitro* Assay to Test sgRNA Cleavage Efficiency

The **Guide-it sgRNA Screening Kit** is a complete system for predicting the cleavage efficacy of sgRNAs *in vitro*, prior to use for genome editing in cells (Figure 1). With this kit, a template containing a sgRNA-target site is created by PCR; then the test sgRNA and recombinant Cas9 nuclease are added. The efficiency of Cas9-mediated cleavage can be measured by agarose gel electrophoresis.


## 1 Use PCR to generate a target for cleavage



## 2 In vitro cleavage of target sequence by recombinant Cas9 and synthesized sgRNA



## 3 Separate cleavage products on an agarose gel




**Figure 1. Overview of the Guide-it sgRNA Screening Kit protocol.** A PCR amplicon containing a sgRNA target site is synthesized from genomic DNA (Step 1). The PCR fragment is then combined with a candidate sgRNA and recombinant Cas9 (Step 2). The entire reaction is separated by agarose gel electrophoresis (Step 3). Since the sgRNA-target sequence is located asymmetrically within the amplicon, cleavage by the Cas9-sgRNA complex results in two bands of unequal length that can be easily distinguished on an agarose gel.

### sgRNAs Exhibit Different Cleavage Efficiencies

CRISPR/Cas9 genome editing was used to disrupt the *CXCR4* locus in HeLa cells. *CXCR4* encodes a cell surface chemokine receptor that interacts with the CXCL12 chemokine and plays an important role in the immune system.


In this experiment, four different sgRNAs targeting the *CXCR4* locus were tested using the [Guide-it sgRNA Screening Kit](#). Briefly, sgRNAs targeting the *CXCR4* gene were synthesized using the [Guide-it sgRNA In Vitro Transcription Kit](#). A PCR fragment containing the sgRNA target sequence was mixed with recombinant Cas9 protein and each sgRNA. The cleavage reaction was analyzed by agarose gel electrophoresis. Densitometry (Cong *et al.*, 2013) showed that sgRNA3 had the lowest cleavage efficiency (Figure 2).



**Figure 2. Differences in *in vitro* cleavage efficiency as determined by the Guide-it sgRNA Screening Kit.** The cleavage efficiency of four different sgRNAs targeting the *CXCR4* locus were tested. A PCR fragment containing the *CXCR4* target sequence was synthesized and mixed with Cas9 and each sgRNA. A negative control that lacked sgRNA was included for comparison (NC). Cleavage efficiency was assessed by agarose gel electrophoresis and measured using densitometry (%).


### *In Vitro* Cleavage Efficiency Predicts *In Vivo* Cleavage

HeLa cells were co-transfected with plasmids encoding Cas9 and one of the four different sgRNAs tested above. The presence of mutations in the *CXCR4* locus as was assayed using the [Guide-it Mutation Detection Kit](#). This assay uses a mismatch-specific nuclease, Guide-it Resolvase, to identify insertions or deletions in specific loci in cells treated with engineered nucleases. Mismatches were detected with high efficiency in cells treated with sgRNA1, 2, and 4 (Figure 3). However, cells treated with sgRNA3 exhibited a very low efficiency of mismatches, consistent with the efficiency predicted by the Guide-it sgRNA Screening Kit (Figure 2).



**Figure 3. sgRNA-mediated cleavage in HeLa cells as determined by the Mutation Detection Kit.** HeLa cells were co-transfected with plasmids encoding Cas9 and one of the four different sgRNAs using [Xfect Transfection Reagent](#). Six days after transfection, cells were assayed for the presence of mutations using Guide-it Resolvase, a mismatch-specific nuclease. Cleavage fragments were present for all sgRNAs except sgRNA3, indicating low Cas9 guiding efficiency for this particular sgRNA.


*CXCR4* gene disruption was also assessed by flow cytometry; since *CXCR4* is a cell surface receptor, it can be detected by flow cytometry using a FITC-labeled *CXCR4* antibody. Disruption in *CXCR4* expression could be detected in cells transfected with Cas9 and sgRNA1, 2, and 4 (Figure 4). In contrast, for cells transfected with Cas9 and sgRNA3, a much smaller proportion of the cells had disruption of *CXCR4* expression. These functional data confirm the results obtained by both the Guide-it sgRNA Screening Kit and Mutation Detection Kit.



**Figure 4. Flow cytometric analysis detects sgRNA-mediated loss of *CXCR4* function.** HeLa cells cotransfected with plasmids encoding Cas9 and one of four sgRNAs, were stained with a FITC-labeled antibody against *CXCR4*. Knockout of *CXCR4* gene by CRISPR/Cas9 editing will result in reduced protein expression. Therefore, FITC staining is inversely correlated with efficient genome editing. The percentage (%) of the cell population that were not labeled with FITC is shown in blue. Cells treated with Cas9 and sgRNA3 exhibit the greatest percentage of FITC+ cells and the least efficient genome editing.

## Conclusions: Accurate Prediction of sgRNA Cleavage Efficiency

There is a clear correlation between *in vitro* sgRNA cleavage efficiency as predicted by the Guide-it sgRNA Screening Kit, and *in vivo* sgRNA-mediated cleavage as assessed by the presence of indels and functional gene knockout (Figure 5). These results indicate that the Guide-it sgRNA Screening Kit is ideal method for screening for ineffective sgRNAs during CRISPR/Cas9 genome editing projects.



**Figure 5. Guide-it sgRNA Screening Kit accurately predicts *in vivo* sgRNA efficacy.** Cleavage efficiency was assessed by *in vitro* cleavage (Figure 2) and the Guide-it Mutation Detection Kit (Figure 3); functional knockout was assessed by flow cytometry (Figure 4, % of CXCR4<sup>-</sup> cells). There is a clear correlation between the efficiency predicted by the Guide-it sgRNA Screening Kit, and estimation of *in vivo* cleavage (Mutation Detection Kit) and functional knockout (Flow Cytometry).

## References

Cong, L. *et al.* (2013) Multiplex genome engineering using CRISPR/Cas9 systems. *Science* **339**(6121):819–23.

[View web page >>](#)

[http://www.clontech.com/US/Products/Genome\\_Editing/CRISPR\\_Cas9/Technical\\_Notes/Screening\\_for\\_effective\\_guideRNAs](http://www.clontech.com/US/Products/Genome_Editing/CRISPR_Cas9/Technical_Notes/Screening_for_effective_guideRNAs)

This is a reprint from a page on our web site. All license, copyright, and trademark information pertaining to this content applies as stated in the original web content. This information can be found at [www.clontech.com](http://www.clontech.com).