Protein Molecular Weight Marker (Broad)

Code No. 3452

Size: for 200 lanes

Conc.: 18 $\mu\text{g}/\mu\text{l}$

Supplied Reagents:

5X Loading Buffer

1 ml

1 M DTT (Dithiothreitol)

100 μl

Form:

50 mM	Tris-HCl, pH 6.8
1 mM	EDTA
200 mM	NaCl
50%	Glycerol

Description:

This marker consists of 9 kinds of proteins with a Molecular Weight range between 6.5 kDa and 200 kDa. Each protein is proportioned to yield uniform band intensities when stained with Coomassie Brilliant Blue R-250 after run on SDS-polyacrylamide gel.

5 μl of 20-fold diluted marker is loaded per lane of SDS-PAGE minigel for standard use. In this case, this marker is sufficient for 200 lanes.

Component Proteins:

Protein	Source	MW (Da)
Myosin	Pig	200,000
β -galactosidase	<i>E. coli</i>	116,000
Phosphorylase B	Rabbit muscle	97,200
Serum Albumin	Bovine	66,409
Ovalbumin	Hen egg White	44,287
Carbonic anhydrase	Bovine	29,000
Trypsin inhibitor	Soybean	20,100
Lysozyme	Hen egg white	14,300
Aprotinin	Bovine pancreas	6,500

Storage:

Protein MW Marker (Broad), 1 M DTT: -20°C

5X Loading Buffer: Store at RT after used

* 2 years from date of receipt under proper storage conditions.

5X Loading Buffer (Store at RT after used):

200 mM	Tris-HCl, pH 6.8
10%	SDS
0.05%	BPB
50%	Glycerol

Application Example:


- Combine the following reagents in a tube.

1 M DTT	2 μl
5X Loading Buffer	20 μl
- Prepare 20-fold diluted marker by adding the following components to the solution prepared in 1.

Protein MW Marker (Broad)	5 μl
Sterile purified water	73 μl

* 20-fold diluted marker is stable for 2-3 months for -20°C. It is recommended to store in aliquots for several uses, not to repeat freeze-thaw cycles.
- Mix 20-fold diluted marker well, and heat at 100°C for 5 minutes. Load 5 μl per lane of SDS-PAGE minigel. Run 5-20% gradient SDS-polyacrylamide gel electrophoresis.
- Perform staining with Coomassie Brilliant Blue R-250.

Electrophoresis Result (5 - 20% gradient SDS-PAGE)

Note:

It is recommended to use 5-20% gradient SDS-polyacrylamide gel for electrophoresis.

If the concentration of gel is too low, the mobility of bands of low molecular weight is faster than the mobility of BPB. If the concentration of gel is too high, high molecular weight bands may not separate well or they may not migrate.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Protein Molecular Weight Marker (Broad)

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