pPTR I DNA

Code No. 3621

Size: 20 μ g

Conc.: 1 μ g/ μ l

* 2 years from date of receipt under proper storage conditions.

Description:

pPTR I is a chromosomal integrating and *E. coli-Aspergillus* shuttle vector. In *Aspergillus* cells, pPTR I does not replicate autonomously and is maintained only when integrated into chromosome. This vector has Ampicillin-resistant gene, Amp^r, as the selection marker for *E. coli* and Pyrithiamine-resistant gene, *ptrA*, as that for *Aspergillus*. Since it has cloning sites, containing *Hind* III, *Pst* I, *Sma* I and *Kpn* I, on *lacZ* region, any DNA insertion into this vector can be easily verified using plates containing IPTG and X-Gal.

Form: 10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage: -20°C

Preparation: Purified by column treatment

Base pairs: 4,782 bp

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Usage:

This is a vector for transformation of *Aspergillus* species such as *A. oryzae*, *A. niger*, *A. terreus* and *A. fumigatus* using Pyrithiamine which is a metabolic competitive analog of thiamine.


Reference:

Kubodera T, Yamashita N, and Nishimura A.
Biosci Biotechnol Biochem. (2000) **64**: 1416-1426.

Note:


As the presence of thiamine causes the remarkable decrement of selection efficiency by Pyrithiamine, the media and the reagents without thiamine should be used for this transformation system.

Vector map of pPTR I:

ptrA : Pyrithiamine resistant gene in *A. oryzae*
lacZ : β -galactosidase gene in *E. coli*
Amp^r : Selective marker in *E. coli*
Coli ori : Replication origin in *E. coli*

Cloning site of pPTR I:

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

This product is sold under licensing with HAKUTSURU SAKE BREWING CO.,LTD.

Pyrithiamine-resistant Transformation Systems

Transformation of *A. oryzae* through Protoplast-PEG method

Sterilize all instruments and solutions to be used, before starting the protocol.

All operations must be performed under sterile condition.

[Preparation]

• Spore suspension :

Collect the spores of *A. oryzae* on a plate and suspend in 10 ml of 0.1% Tween 80 and 0.8% NaCl. Filter the suspension with a glass filtration (3G2) to remove the hypha and the agar debris, and recover the filtered spore solution.

Centrifuge at 3,000 rpm for 5 min to collect the spores, and discard the supernatant. Wash the spores with 10 ml of 0.1% Tween 80 twice. Suspend the spores in a proper amount of sterile purified water.

• CD medium (per 1 liter) :

| | |
|---|-----------|
| NaNO ₃ | 6.0 g |
| KCl | 0.52 g |
| KH ₂ PO ₄ | 1.52 g |
| Glucose | 10.0 g *1 |
| 1 M MgSO ₄ · 7H ₂ O | 2 ml *2 |
| Trace elements solution *3 | 1 ml |
| Agar | 20.0 g) |

Adjust to pH 6.5 with 1 N KOH.

* 1 MgSO₄ · 7H₂O should be autoclaved separately and mixed later.

* 2 When a strain shows very bad growth, culture at the higher glucose concentration around up to 20.0 gram.

* 3 Trace elements solution (per 1 liter)

| | |
|---|--------|
| FeSO ₄ · 7H ₂ O | 1.0 g |
| ZnSO ₄ · 7H ₂ O | 8.8 g |
| CuSO ₄ · 5H ₂ O | 0.4 g |
| Na ₂ B ₄ O ₇ · 10H ₂ O | 0.1 g |
| (NH ₄) ₆ Mo ₇ O ₂₄ · 4H ₂ O | 0.05 g |

• CD selective medium :

CD medium added with 0.8 M NaCl and 0.1 μg/ml Pyrithiamine (PT).

• CD soft agar selective medium :

CD selective medium having 0.5% agar.
Keep it at approximately 50°C.

• Protoplast forming solution : Sterilize by filtering.

20 mg/ml Yatalase (Cat. #T017), 0.8 M NaCl,
10 mM Na phosphate buffer (pH 6.0).

• Protoplast washing solution : 0.8 M NaCl

• Solution 1 :

0.8 M NaCl, 10 mM CaCl₂, 10 mM Tris-HCl (pH 8.0)

• Solution 2 : Sterilize by filtering.

40% (w/v) PEG4000, 50 mM CaCl₂, 50 mM Tris-HCl (pH 8.0)

[Protocol]

- 1) Inoculate the spore suspension of *A. oryzae* in 100 ml of CD liquid medium. *4 Culture by shaking at 30°C for 20 hours. *5 (Shake vigorously in a triangle flask with baffle.)
 - 2) Filter with a glass filtration (3G1) to collect the hypha. Wash with sterile purified water and remove water from the hypha fully by pressing with a spatula.
 - 3) Suspend proper amount of hypha in protoplast forming solution contained in a 50 ml polypropylene centrifuge tube.
 - 4) Shake gently at 30°C and form protoplast. Confirm the formation of protoplast by observing the status with a microscope every 30 min. (Usually the reaction is completed within 2 hours.)
 - 5) Filter with a glass filtration (3G2). Centrifuge the filtered solution at 2,000 rpm for 5 min to collect protoplast.
 - 6) Wash the protoplast twice with 0.8 M NaCl.
 - 7) Suspend the protoplast in Solution 1 to have the concentration of 2×10^8 /ml. Then add the 20% volume of Solution 2 and suspend gently. As Solution 2 has high viscosity, suspend the solution fully.
 - 8) Add plasmid (in less than 20 μl) into 0.2 ml of protoplast suspension contained in polypropylene centrifuge tube (larger than 10 ml). Plasmid can be added at maximum 20 μg.
 - 9) Leave to stand on ice for 30 min.
 - 10) Add 1 ml of Solution 2. Suspend gently.
 - 11) Leave to stand at room temperature for 15 min.
 - 12) Add 8.5 ml of Solution 1. Suspend gently.
 - 13) Centrifuge the solution to collect the protoplast. Discard the supernatant as much as possible. Suspend the collected protoplast in 0.2 ml of Solution 1.
 - 14) Mix the protoplast suspension in 5 ml of CD soft agar selective medium. *4 Quickly overlay the protoplast/soft agar on CD selective plate (90 mm-diameter) *4 to be spread uniformly.
 - 15) Culture at 30°C for 5 - 7 days.
- * 4 Use CD medium without thiamine. Under the existence of thiamine, the selection efficiency by PT drastically decreases.
- * 5 Some strains may allow good results after cultured for 40 - 48 hours by shaking.

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