Adenovirus genome DNA-TPC

Code No. 6171

Size : 50 μ l (5 tests)

Shipping at – 80°C

Store at – 80°C

Note:

The efficiency may be decreased by freeze-thawing of DNA-TPC.
Avoid repeated freeze-thaw cycles unnecessarily.

Lot No.

Volume : 50 μ l (5 tests)

Expiry Date :

Note:

For construction of recombinant adenovirus, both Adenovirus Adenovirus Dual Expression Kit (Cat. #6170) and this product are required.

Description :

The COS-TPC method¹⁾ is a highly efficient adenovirus construction method developed by Dr. I. Saito, et al., the University of Tokyo. This product is used to construct recombinant adenovirus through the COS-TPC method by being transfected in combination with the recombinant cosmid which is prepared with Takara's Adenovirus Adenovirus Dual Expression Kit (Cat. #6170). This product is attached with TP (Terminal Protein) which is originally bound to both termini of adenovirus genome DNA. Accordingly, the COS-TPC method achieves the highly efficient construction of recombinant adenovirus, and the target clones are produced with higher efficiency. This product was already digested with *Eco*T22 I and *Cla* I.

Related Products :

Adenovirus Adenovirus Dual Expression Kit (Cat. #6170)

Protocol :

This product is used by following the protocol which is supplied with Adenovirus Adenovirus Dual Expression Kit (Cat. #6170).

Reference :

1) Miyake,S, Makimura,M., Kanegae,Y., Harada,S., Takamori,K., Tokuda,C., and Saito,I. (1996) *Proc. Natl. Acad. Sci. USA*, **93**, 1320

Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at www.takara-bio.com.

The viral supernatants produced by this retroviral vector could, contain potentially hazardous recombinant virus depending on an insert gene. Due caution must be exercised in this production and handling of recombinant retrovirus. Please follow the guideline for experiments using recombinant DNA issued by the relevant authorities and the safety committee of your organization or your country in using this product.

Takara is not liable for any accident or damage caused by the use of this product.

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