# Lysis Buffer for *Legionella*

Code No. 9181

Size: 50 preps

## Component:

Lysis Buffer for *Legionella* 500  $\mu$ l x 5

### Description:

This product is a DNA extraction reagent for *Legionella* species, composed with detergent and PCR inhibitor-absorbing resin. DNA is extracted by heat treatment and centrifugation. As the gel layer is formed on the resin after centrifugation, DNA solution can be collected easily.

**Storage:** 4°C (do not freeze)


### Protocol:

1. Centrifuge the solution containing *Legionella* bacteria for 5 minutes at 15,000 rpm (maximum speed) at 4°C. Remove and discard the supernatant.
2. Mix the Lysis Buffer for *Legionella* well by vortexing before use. Pipette to mix prior to addition to the sample to ensure uniform mixing of the resin. Add 50  $\mu$ l of Lysis Buffer for *Legionella* to the cell pellet. Mix briefly by vortexing and spin down briefly.
3. Incubate at 95°C for 10 minutes.
4. Mix briefly by vortexing. Centrifuge the solution at 15,000 rpm (maximum speed) at 4°C for 10 minutes.
5. Incubate on ice for 5 minutes.
6. Recover the supernatant. Approximately 25  $\mu$ l of DNA solution should be recovered.

\* A low concentration gel has been added to the Lysis Buffer for *Legionella*. After cell lysis and centrifugation, a gel layer is formed on the resin (the gel layer can not be seen visually). In step 5., incubation on ice firmly solidifies the gel layer. To ensure that the tip does not touch the gel layer when collecting the supernatant, pipette from near the DNA solution surface.

### Notes:

- Mix well by vortexing before use. Also, mix by pipetting prior to fractionation to ensure uniform distribution of the resin. In addition, if tip clogging occurs, cut the tip to enlarge the opening.

### ◆ Note ◆

Do not touch the tip to the gel layer; collect DNA solution from near the surface of the supernatant.

### Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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