# pLVpro2-EF1 $\alpha$ -Km Vector

Code No. 6979

Size: 20  $\mu$ g

Conc.: 0.5  $\mu$ g/ $\mu$ l

## Description:

pLVpro2-EF1  $\alpha$ -Km Vector is a pLVpro2 lentiviral vector plasmid that can be used by inserting a target gene sequence (cDNA or gene fragment) into a Multi Cloning Site (MCS).

pLVpro2 Lentivirus Vector Plasmid is a SIN-type lentiviral vector plasmid that contains the sequences that, together with the HIV-1 LTR (CMV- 5'LTR and 3'LTR/ $\Delta$  U3) and the Lentivirus packaging signal ( $\Psi$ ), can improve the expression of the target gene, the virus titer, and the overall vector function.

By simply co-transfect Lenti-X™ 293T cells (Lenti-X 293T Cell Line, Cat. #632180) with pLVpro2 Lentivirus Vector Plasmid carrying your gene of interest, and LVpro™ Packaging Mix (Cat. #6195), non-replicating recombinant virions (virus particles) can be easily prepared.

\* For more information on the LVpro Packaging Mix, see the user manual for LVpro Packaging Mix with pLVpro Series (Cat. #6195, 6962 - 6967).

Storage: -20°C

Base Pairs: 6,823 bp


## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## MCS:

Sse8387 I      Not I      Sml I  
BamH I      Hpa I      Mlu I      Xba I  
GGATCCTGCAGGTTAACCGCGCCGCACGCGTCTCGAGGGCCCATCGATTTAAATCTAGA  
CCTAGGACGTCCAATTGCCCGCGCTGCGCAGAGCTCCCGGTAGCTAAATTTAGATCT

## Vector Map of pLVpro2-EF1 $\alpha$ -Km Vector:

## Biosafety:

The protocols in this User Manual require the production, handling, and storage of infectious lentivirus. It is imperative to fully understand the potential hazards of, and take necessary precautions for the laboratory use of lentiviruses.

Although biosafety regulations and practices vary from country to country, the National Institute of Health and Center for Disease Control and Prevention (CDC) have designated recombinant lentiviruses as Level 2 organisms. This requires the maintenance of a Biosafety Level 2 facility for work involving this virus and others like it.

The pseudotyped lentiviruses packaged from the HIV-1-based vectors described here are capable of infecting human cells. The viral supernatants produced by these lentiviral systems could, depending on your insert, contain potentially hazardous recombinant virus. Similar vectors have been approved for human gene therapy trials, attesting to their potential ability to express genes *in vivo*.

For more information on Biosafety Level 2 agents and practices, download the following reference: Biosafety in Microbiological and Biomedical Laboratories (BMBL), 6th Edition (2020) HHS Pub. No. (CDC) 300859. U.S. Department of Health and Human Services CDC and NIH. If possible, observe and learn the practices described below from someone who has experience working with lentiviruses.

Summary of Biosafety Level 2: Practices: - Standard microbiological practices - Limited access to work area - Biohazard warning signs posted - Minimize production of aerosols - Decontaminate potentially infectious wastes before disposal - Use precautions with sharps (e.g., syringes, blades) - Biosafety manual defining any needed waste decontamination or medical surveillance policies.

**Note:** Please be advised that our company cannot be held responsible for any accidents or damage related to the handling and use of this product.

LVpro is a trademark of Takara Bio Inc.

Lenti-X is a trademark of Takara Bio USA, Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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