Legionella LC Medium Base

Code No. 9016 Size: 90 ml

Description:

This product is a base of the medium for liquid culture (LC) in LC EMA-qPCR method for *Legionella* species. It is used by adding the commercially available supplements for *Legionella*.

Composition (in 90 ml):

Yeast Extract	1.0 g
Activated charcoal (Norit SA2)	0.2 g

Storage: -20°C

Materials required but not provided:

Legionella Growth Supplement (BCYE)
(KANTO CHEMICALS (Oxoid), Code. 713251-1)
Legionella Selective Supplement (MWY)
(KANTO CHEMICALS (Oxoid), Code. 713255-1)

[NOTE] Please contact us for the preparation of MWY media when using the supplement of the previous version (KANTO CHEMICALS (Oxoid); Code. 713110-1 and 713118-1).

Preparation of MWY Media:

1. Bring *Legionella* LC Medium Base to room temperature.
2. Add 0.05 g of α -ketoglutarate and 1.0 g of ACES buffer to 90 ml of *Legionella* LC Medium Base and dissolve.
3. Add 3.52 ml of 1 M KOH. *¹
4. Dissolve *Legionella* Growth Supplement (BCYE) in 5 ml of sterile purified water (Hot water: less than 50°C). Then add all 5 ml to the solution from step 3, and mix.
5. Dissolve *Legionella* Selective Supplement (MWY) in 2 ml of sterile purified water. Then add all 2 ml to the solution from step 4, and mix.
6. Aliquot 1 ml in sterile microtubes and store at -20°C. *²

*¹ At this point, the pH of the solution becomes 7.1 ± 0.1 . After the supplement is added, the pH of the solution finally becomes 6.9 ± 0.1 .

*² Microtip clogging with activated charcoal may occur.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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