Anti-Human Influenza A (H1N1, H2N2), Monoclonal (Clone C179)

Code No. M145

Size : 0.1 mg Mouse Ig

Subclass : IgG2a

* 2 years from date of receipt under proper storage conditions.

Source :

This monoclonal antibody was obtained by fusing the mouse myeloma cell line P3U1 with spleen cells of a BALB/c mouse immunized with influenza A virus A/Okuda/57 strain (H2N2 subtype). The monoclonal antibody was harvested from ascites fluid.

Purification :

The antibody was purified by column chromatography, dissolved in 10 mM PBS pH 7.4 containing 1.0% bovine serum albumin, and then lyophilized. The lyophilized antibody does not contain preservatives.

Form : Lyophilized

Reconstitution :

Dissolve the lyophilized antibody in 50 μ l of distilled water (final antibody concentration: 2.0 mg/ml). This solution can be used as a stock solution. If dilution is necessary, dilute the stock solution with the Dilution solution described below immediately prior to use. When the entire amount of antibody will be used over a short time period, it may be dissolved directly in 500 μ l or more of the Dilution solution.

Note (1): Store the antibody at a minimum concentration of 2.0 mg/ml. Lower antibody concentration may decrease stability.

Note (2): The reconstituted antibody solution should contain 0.1% sodium azide as a preservative when stored at 4°C.

Dilution solution: 10 mM PBS (pH 7.4)
1.0% BSA
(0.1% NaN₃)*

* When stored at 4°C, 0.1% sodium azide should be added as a preservative.

Specificity:

- This antibody specifically recognizes the conformational epitope in the stem region of HA; the amino acid sequence of this region is highly conserved in H1N1 and H2N2 subtypes of influenza A virus.
- This antibody does not cross-react with HA in the H3N2 subtype of influenza A virus.
- Refer to the Table 1 for reactivity of monoclonal antibodies against influenza A and B viruses.
- This antibody has neutralization activity for H1N1 and H2N2 subtypes.
- This antibody reacts with the H5N3 subtype of chicken influenza.
- This antibody has no effect on hemagglutination-inhibition for all subtypes of virus.

Storage : 4°C

This product does not contain a preservative.

Store the stock solution (2.0 mg/ml) in aliquots at -20°C for 1 year, or at 4°C for 6 months after adding 0.1% sodium azide. Avoid repeated freezing and thawing. Diluted antibody should not be stored.

Working concentration :

- 0.5 - 1 μ g/ml (For immunohistochemical staining by the PAP [peroxidase-antiperoxidase] staining method)
- 10 μ g/ml (For neutralization test)

Application :

- Influenza virus typing
- Neutralization activity test of influenza virus H1N1 and H2N2
- Protection against influenza virus in animal infection experiments

Caution: This antibody can not be used for Western Blotting.

References :

- 1) Okuno Y. *J Virol.* (1993) **67**(5): 2552-2558.
- 2) Ueda M and Okuno Y. *J Clin Microbiol.* (1998) **36**(2): 340-344.
- 3) Okuno Y. *J Virol.* (1994) **68**(1): 517-520. De nobit aute del il ipsaevolorepre prae modit ut fugit di volo que prempore volum assitas aces

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Anti-Human Influenza A (H1N1, H2N2), Monoclonal (Clone C179)

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