

**pMEI-5 DNA(製品コード 3656)**

Vector Elements	Location	由来する生物種	実験分類 クラス
5' LTR MLV U3	408...853	Murine leukemia virus	2
5' LTR MLV R	854...922	Murine leukemia virus	2
5' LTR MLV U5	923...999	Murine leukemia virus	2
SD	1057...1062	Murine leukemia virus	2
Ψ (packaging signal)	1069...1419	Murine leukemia virus	2
Intron+SA (EF1α intron)	1482...1715	Homo sapiens	1
MCS (multi cloning site)	1760...1816	人工合成核酸	—
3' LTR MLV U3	1861...2306	Murine leukemia virus	2
3' LTR MLV R	2307...2375	Murine leukemia virus	2
3' LTR MLV U5	2376...2452	Murine leukemia virus	2
pUC ori (pUC origin of replication)	2827...3481	<i>Escherichia coli</i>	1
Ampr (Ampicillin resistance gene)	3629...4489	<i>Escherichia coli</i>	1