

## Tissue Microarray – Normal adult and foetal bone tissue

Cat.-No.: 401 1211 Lot: 069P

	a	b	c	d
1	•	•	•	•
2	•	•		
3	•	•		
4	•	•	•	•
5	•	•	•	•
6	•	•	•	•

Technical Information: 20 Spots

- Spot diameter: 2.0 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Localization	Tissue	Sex	Age
1a	radius, ulna	fetal		18.SSW
1b	radius, ulna	fetal		18.SSW
1c	radius, ulna	fetal		18.SSW
1d	radius, ulna	fetal		18.SSW
2a	femur	fetal		18.SSW
2b	femur	fetal		18.SSW
3a	radius, ulna	fetal		18.SSW
3b	radius, ulna	fetal		18.SSW
4a	hip	adult	W	63
4b	hip	adult	W	65
4c	pallus, clacaneus	adult	W	65
4d	pallus, clacaneus	adult	W	65
5a	hip	adult	W	92
5b	hip	adult	W	92
5c	hip	adult	W	62
5d	hip	adult	W	62
6a	hip	adult	M	37
6b	hip	adult	M	37
6c	hip	adult	W	90
6d	hip	adult	W	90

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.