

Tissue Microarray - Cartilage I+II

Human Normal Adult and Neonatal Cartilage Tissue

Cat.-No.: 401 1223 Lot: 154P230910

	a	b	c	d	e
1	•	•	•	•	
2	•	•	•	•	
3	•	•	•	•	
4	•	•	•	•	
5	•	•	•	•	
6	•	•	•	•	
7	•	•	•	•	
8	•	•	•	•	
9	•	•	•	•	
10	•	•			

Technical Information: 38 spots

- Spot diameter: 2.0 mm
- Fixation in 4 % paraformaldehyde in PBS
- Paraffin embedded
- Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	E-Number	Localisation	tissue	Sex	Age
1a	S209-03 K172	hip	adult	male	73
1b	S209-03 K172	hip	adult	male	73
1c	S197-03 K171	hip	adult	male	79
1d	S197-03 K171	hip	adult	male	79
2a	S193-03 K170	hip	adult	female	80
2b	S193-03 K170	hip	adult	female	80
2c	S192-03 K169	hip	adult	female	50
2d	S192-03 K169	hip	adult	female	50
3a	A21-02 K10	knee	adult	male	76
3b	A21-02 K10	knee	adult	male	76
3c	A39-02 K12	knee	adult	male	56
3d	A39-02 K12	knee	adult	male	56
4a	E29291-07	knee	adult	female	45
4b	E29291-07	knee	adult	female	45
4c	E12369-07	knee	adult	female	33
4d	E12369-07	knee	adult	female	33
5a	A116-02 K21	knee	adult	female	83
5b	A116-02 K21	knee	adult	female	83
5c	A137-02 K20	Knee (synovia)	adult	female	52
5d	A137-02 K20	Knee (synovial)	adult	female	52
6a	S698-01 K25b	humerus	fetal, normal		30 weeks of gestation
6b	S698-01 K25b	humerus	fetal, normal		30 weeks of gestation
6c	S698-01 K25a	femur	fetal, normal		30 weeks of gestation
6d	S698-01 K25a	femur	fetal, normal		30 weeks of gestation
7a	S739-00 K28	femur re	fetal, normal		23 weeks of gestation

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	E-Nummer	Localisation	tissue	Sex	Age
7b	S739-00 K28	femur re	fetal, normal		23 weeks of gestation
7c	S960-00 K29	femur + tibia	fetal, normal		27 weeks of gestation
7d	S960-00 K29	femur + tibia	fetal, normal		27 weeks of gestation
8a	S481-02 K30	femur prox.	fetal, normal		27 weeks of gestation
8b	S481-02 K30	femur prox.	fetal, normal		27 weeks of gestation
8c	E33281-02 K119 II	humerus	fetal, normal		13 weeks of gestation
8d	E33281-02 K119 II	humerus	fetal, normal		13 weeks of gestation
9a	E20034-02 K95b		embryonal, normal		12 weeks of gestation
9b	E20034-02 K95b		embryonal, normal		12 weeks of gestation
9c	E7257-02 K47b		embryonal, normal / Plazenta		9 weeks of gestation
9d	E7257-02 K47b		embryonal, normal		9 weeks of gestation
10a	E6701-02 K48d		embryonal, normal		10 weeks of gestation
10b	E6701-02 K48d		embryonal, normal		10 weeks of gestation

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.