

Tissue Microarray - Normal tissue, multi-species

Cat.-No.: 401 1310 Lot: 063P

	Label			
	human a	pig b	rat c	mouse d
1 heart	•	•	•	•
2 lung	•	•	•	•
3 liver	•	•	•	•
4 kidney	•	•	•	•
5 stomach	•	•	•	•
6 small intestine	•	•	•	•
7 colon	•	•	•	•
8 muscle	•	•	•	•
9 skin	•	•	•	•
10 spleen	•	•	•	•
11 brain	•	•	•	•
12 fat	•	•	•	•

Technical Information: 48 Spots

- Spot diameter: 1.5 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded
- Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Species	Breed	Organ	Sex
1a	human		heart	F
1b	domestic pig		heart	M
1c	inbred rat	lewis	heart	M
1d	Immunodeficient mouse	CD-1® Nude	heart	F
2a	human		lung	M
2b	domestic pig		lung	M
2c	inbred rat	lewis	lung	M
2d	Immunodeficient mouse	CD-1® Nude	lung	F
3a	human		liver	M
3b	domestic pig		liver	M
3c	inbred rat	lewis	liver	M
3d	Immunodeficient mouse	CD-1® Nude	liver	F
4a	human		kidney	F
4b	domestic pig		kidney	M
4c	inbred rat	lewis	kidney	M
4d	Immunodeficient mouse	CD-1® Nude	kidney	F
5a	human		stomach	M
5b	domestic pig		stomach	M
5c	inbred rat	lewis	stomach	M
5d	Immunodeficient mouse	CD-1® Nude	stomach	F
6a	human		small intestine	M
6b	domestic pig		small intestine	M
6c	inbred rat	lewis	small intestine	M
6d	Immunodeficient mouse	CD-1® Nude	small intestine	F

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Species	Breed	Organ	Sex
7a	human		colon	M
7b	domestic pig		colon	M
7c	inbred rat	lewis	colon	M
7d	Immunodeficient mouse	CD-1® Nude	colon	F
8a	human		muscle	M
8b	domestic pig		muscle	M
8c	inbred rat	lewis	muscle	M
8d	Immunodeficient mouse	CD-1® Nude	muscle	F
9a	human		skin	F
9b	domestic pig		skin	M
9c	inbred rat	lewis	skin	M
9d	Immunodeficient mouse	CD-1® Nude	skin	F
10a	human		spleen	F
10b	domestic pig		spleen	M
10c	inbred rat	lewis	spleen	M
10d	Immunodeficient mouse	CD-1® Nude	spleen	F
11a	human		brain	M
11b	domestic pig		brain	M
11c	inbred rat	lewis	brain	M
11d	Immunodeficient mouse	CD-1® Nude	brain	F
12a	human		fat	F
12b	domestic pig		fat	M
12c	inbred rat	lewis	fat	M
12d	Immunodeficient mouse	CD-1® Nude	fat	F

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.