

Tissue Microarray – Stem cell rich tissue - Stem Cell TMA[®]

Cat.-No.:401 1401 Lot: 054P

	a	b	c	d
1	•	•	•	
2	•	•		
3	•	•	•	•
4	•	•	•	
5	•	•	•	
6	•	•	•	

Technical Information: 18 Spots

- Spot diameter: 2.0 mm
 - Fixation in 4% paraformaldehyde in PBS
 - Paraffin embedded
- Tissue type validated by immunohistochemistry CD 34 +

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Tissue	Sex	Week of gestation
1a	abdomen	n.s.	12.
1b	abdomen	n.s.	12.
1c	liver	m	34.
2a	liver	m	34.
2b	liver	m	21.
3a	liver	m	21.
3b	liver	f	25.
3c	liver	f	25.
3d	liver	f	32.
4a	liver	f	32.
4b	liver	m	17.
4c	liver	m	17.
5a	liver	f	18.
5b	liver	f	18.
5c	liver	n.s.	14.
6a	liver	n.s.	14.
6b	liver	n.s.	14.
6c	liver	n.s.	14.

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.