

Tissue Microarray - Lymphoma

Cat.-No.: 401 2101 Lot: 018P

	a	b	c	d	e
1	•	•			
2	•	•	•		
3	•	•	•	•	•
4	•	•	•	•	•
5	•	•	•	•	•

Technical Information: 20 spots

- Spot diameter: 2mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Tissue	Sex	Age	Diagnostics	ICD-O
1a	Lymph node	male	53	normal	
1b	Lymph node	male	53	normal	
2a	tonsil	male	51	normal	
2b	tonsil	male	51	normal	
2c	Lymph node	male	64	BCLL	9823/3
3a	Lymph node	male	64	BCLL	9823/3
3b	Lymph node	male	63	Mantle cell lymphoma	9673/3
3c	Lymph node	male	63	Mantle cell lymphoma	9673/3
3d	Lymph node	male	60	DLBCL	9680/3
3e	Lymph node	male	60	DLBCL	9680/3
4a	Lymph node	male	39	DLBCL	9680/3
4b	Lymph node	male	39	DLBCL	9680/3
4c	Lymph node	female	69	DLBCL	9680/3
4d	Lymph node	female	69	DLBCL	9680/3
4e	Lymph node	female	33	Hodgkin	9652/3
5a	Lymph node	female	33	Hodgkin	9652/3
5b	Lymph node	female	20	Hodgkin	9663/3
5c	Lymph node	female	20	Hodgkin	9663/3
5d	Lymph node	male	62	T-NHL-NOS	9702/10
5e	Lymph node	male	62	T-NHL-NOS	9702/10

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.