

Tissue Microarray - Human Mamma Carcinoma

Cat.-No.: 401 2202 Lot: 016P

	a	b	c	d	e	f
1	•	•	•	•	•	•
2	•	•	•	•	•	•
3	•	•	•	•	•	•
4	•	•	•	•	•	•
5	•	•	•	•	•	•
6	•	•	•	•	•	•
7	•	•	•	•		
8	•	•	•	•		•
9	•	•	•	•	•	•
10	•	•	•			
11	•	•	•	•	•	

Technical Information: 59 spots

- Spot diameter: 1,5 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Organ (Specifications)	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
1a	ID-CA	female	67	2	0(0/12)	x	2
1b	ID-CA	female	90	1c	1biii	x	2
1c	ID-CA	female	25	2	x(0/2)	x	3
1d	ID-CA	female	49	2	1bi(1/18)	x	1
1e	ID-CA	female	60	3	1biii	x	3
1f	ID-CA	female	56	1c	0	x	1
2a	ID-CA	female	58	2	x	x	3
2b	ID-CA	female	39	1c	0	x	1
2c	ID-CA	female	73	2	0	x	3
2d	ID-CA	female	56	2	1biv	x	2
2e	ID-CA	female	42	2	1biii	x	3
2f	ID-CA	female	42	1c	1biv	x	3
3a	ID-CA	female	52	1c	0	x	2
3b	ID-CA	female	54	1c	0	x	1
3c	ID-CA	female	66	1c	1biii	x	1
3d	ID-CA	female	53	1b	0	x	2
3e	ID-CA	female	72	2	1biii	x	1
3f	ID-CA	female	46	1c	x	x	1
4a	ID-CA	female	52	1b	x	x	3
4b	ID-CA	female	40	2	x	x	2
4c	ID-CA	female	58	2	0	x	2
4d	ID-CA	female	60	4b	0	x	3
4e	ID-CA	female	63	1c	1biv(9/54)	x	2
4f	ID-CA	female	62	2	1a	x	2
5a	ID-CA	female	69	1b	0	x	1

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Organ (Specifications)	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
5b	ID-CA	female	70	2	0	x	3
5c	ID-CA	female	45	1c	1biii	x	2
5d	ID-CA	female	48	2	1biv	x	2
5e	ID-CA	female	46	1b	0(0/20)	x	3
5f	ID-CA	female	65	4b	2	x	2
6a	ID-CA	female	65	1c	0	x	3
6b	ID-CA	female	57	1c	0	x	3
6c	ID-CA	female	75	2	x	X	2
6d	ID-CA	female	65	2	0	x	2
6e	ID-CA	female	59	1c	0	x	2
6f	ID-CA	female	50	2	0	x	2
7a	ID-CA	female	64	1c	0	x	1
7b	ID-CA	female	63	1b	0	x	2
7c	ID-CA	female	43	2	2	x	3
7d	ID-CA	female	25	2	0	x	3
8a	IL-CA	female	52	2	0	x	2
8b	IL-CA	female	50	2	1bi	x	2
8c	IL-CA	female	59	2	1biv	x	2
8d	IL-CA	female	56	1c	0	x	2
8e	Missing spot						
8f	IL-CA	female	60	1c	1a	x	1
9a	IL-CA	female	57	1c	0	x	2
9b	IL-CA	female	45	1c	0(0/11)	x	2
9c	IL-CA	female	47	2	1biii(4/15)	x	2
9d	IL-CA	female	55	1b	0(0/15)	x	2

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Organ (Specifications)	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
9e	IL-CA	female	69	1c	0	x	2
9f	IL-CA	female	59	1c	x	x	2
10a	IL-CA	female	65	1c	x	x	2
10b	IL-CA	female	58	2	1a	x	2
10c	IL-CA	female	74	2	0	x	2
11a	normal breast	female	48				
11b	normal breast	female	30				
11c	normal breast	female	39				
11d	normal breast, fibrous tissue	female	21				
11e	normal breast, fibrous tissue	female	61				

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.