

Tissue Microarray - Skeletal Tumors

Cat.-No.: 401 2212 Lot: 070P

	a	b	c	d	e	f
1	•	•	•	•	•	•
2	•	•	•	•	•	•
3	•	•	•	•	•	•
4	•	•	•	•	•	•
5	•	•	•	•		
6	•	•	•	•	•	•
7	•	•	•	•	•	•
8	•	•	•	•	•	•
9	•	•	•	•		
10	•	•	•	•		

Technical Information: 54 Spots

- Spot diameter: 2.0 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Localization	Tissue	Sex	Age
1a	Fibrosarkom	thorax	adult	m	51
1b	Fibrosarkom	thorax	adult	m	51
1c	Fibrosarkom	skin	adult	m	65
1d	Fibrosarkom	skin	adult	m	65
1e	Fibrosarkom	skin	adult	m	65
1f	Fibrosarkom	skin	adult	m	65
2a	Osteosarkom	pelvis	adult	f	47
2b	Osteosarkom	pelvis	adult	f	47
2c	Osteosarkom	maxilla	adult	f	28
2d	Osteosarkom	maxilla	adult	f	28
2e	Chondrosarkom	akromion	adult	m	54
2f	Chondrosarkom	akromion	adult	m	54
3a	Chondrosarkom	maxilla	adult	m	25
3b	Chondrosarkom	maxilla	adult	m	25
3c	Chondrosarkom	pelvis	adult	m	44
3d	Chondrosarkom	pelvis	adult	m	44
3e	solitary bone cyst	tibia	adolescent	f	14
3f	solitary bone cyst	tibia	adolescent	f	14
4a	aneurysmatic bone cyst	tibia	adult	m	38
4b	aneurysmatic bone cyst	tibia	adult	m	38
4c	aneurysmatic bone cyst	tibia	adult	m	38
4d	aneurysmatic bone cyst	tibia	adult	m	38
4e	Bone tissue - border area of aneurysmatic bone cyst	pelvis	adolescent	f	13
4f	Bone tissue - border area of aneurysmatic bone cyst	processus coracoideus	adult	f	38

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Localization	Tissue	Sex	Age
5a	fibr. Dysplasia	maxilla	adult	m	60
5b	fibr. Dysplasia	maxilla	adult	m	60
5c	Enchondrom	humerus	adult	f	48
5d	Enchondrom	humerus	adult	f	48
6a	perisotal Chondrom	os metacarpale I	adult	f	65
6b	perisotal Chondrom	os metacarpale I	adult	f	65
6c	Chondrom	os metacarpale I	adult	f	81
6d	Chondrom	os metacarpale I	adult	f	81
6e	Chondrom	femur	adult	f	54
6f	Chondrom	femur	adult	f	54
7a	Chondrom	os metacarpale I	adult	f	81
7b	Chondrom	os metacarpale I	adult	f	81
7c	Osteom	ethmoid	adult	m	79
7d	Osteom	ethmoid	adult	m	79
7e	Osteom	mandible	adolescent	f	13
7f	Osteom	mandible	adolescent	f	13
8a	Osteom	ethmoid	adult	f	19
8b	Osteom	ethmoid	adult	f	19
8c	Osteom	maxillary sinus	adult	m	40
8d	Osteom	maxillary sinus	adult	m	40
8e	Osteom	anconoid	adult	m	53
8f	Osteom	anconoid	adult	m	53
9a	Osteochondrom	os metacarpale V	adult	m	35
9b	Osteochondrom	os metacarpale V	adult	m	35
9c	Osteochondrom	tibia	adult	f	22
9d	Osteochondrom	tibia	adult	f	22

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Localization	Tissue	Sex	Age
10a	Osteochondrom	shoulder	adult	f	22
10b	Osteochondrom	shoulder	adult	f	22
10c	Osteochondrom	fibula	adolescent	f	17
10d	Osteochondrom	fibula	adolescent	f	17

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.