

Human Tissue Microarray - Multi Tumor : 4 organs

Cat.-No.: 401 2401 Lot: 283P.1.050914.01-05

Slide label				
	a	b	c	d
1	•	•	•	•
2	•	•	•	•
3	•	•		•
4	•	•		•

Technical Information: 12 + 2 additional spots = 14 spots

- Spot diameter: 2 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Organ	Sex	Age	Stage pT	Stage pN	Stage pM	Grade	Remarks
1a	ADC	lung	m	63	2	0	X	2	
1b	ADC	lung	m	66	3	1	X	3	
1c	Peritumorous fibromuscular tissue	colon	f	78	2	0	X	3	mainly necrosis, replacement spot Pos. 4d
1d	ADC	colon	m	68	3	2	X	2	
2a	ADC	prostate	m	73	3a	0	X	3a	
2b	ADC	prostate	m	66	2b	0	X	2a	
2c	ADC	breast	f	65	2	3a	X	2	
2d	ADC	breast	f	65	y4b	y1a	1	2	
3a	normal	lung	f	47					
3b	normal	colon	m	68					Same patient as sample 1d limited spot area,
4a	Normal (hyperplasia)	prostate	m	62					
4b	normal	breast	f	31					
3d	normal	breast	f	31					additional spot
4d	ADC	colon	f	91	4a	1b	x	3	additional spot

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.