

Tissue Microarray - Human Synovitis Tissue

Cat.-No.: 401 3201 Lot: 061P

	a	b	c	d	e
1	•	•		•	•
2	•	•		•	•
3	•	•		•	•
4	•	•	•	•	
5	•	•		•	•
6	•	•		•	•
7	•	•		•	•

Technical Information: 28 Spots

- Spot diameter: 2.0 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded
- Tissue type validated by immunohistochemistry /classified according to synovitis score by Krenn (Pathol Res Pract. 2002 and Pathol Res Pract. 2004; 200(2):95-103)

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Tissue	Diagnosis	Score	Sex	Age
1 a	knee	RA	3 / 9	f	58
1 b	knee	RA	3 / 9	f	58
1 d	knee	RA	4 / 9	f	65
1 e	knee	RA	4 / 9	f	65
2 a	knee	RA	6 / 9	f	57
2 b	knee	RA	6 / 9	f	57
2 d	knee	RA	7 / 9	f	46
2 e	knee	RA	7 / 9	f	46
3 a	knee	RA	8 / 9	f	39
3 b	knee	RA	8 / 9	f	39
3 d	hip	PSA	2 / 9	m	55
3 e	hip	PSA	2 / 9	m	55
4 a	hip	PSA	6 / 9	m	28
4 b	hip	PSA	6 / 9	m	28
4 c	knee	normal tissue	0 / 9	m	24
4 d	knee	normal tissue	0 / 9	m	24
5 a	knee	normal tissue	0 / 9	f	45
5 b	knee	normal tissue	0 / 9	f	45
5 d	knee	OA	0 / 9	f	65
5 e	knee	OA	0 / 9	f	65
6 a	hip	OA	1 / 9	m	74
6 b	hip	OA	1 / 9	m	74
6 d	hip	OA	2 / 9	m	54
6 e	hip	OA	2 / 9	m	54

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Tissue	Diagnosis	Score	Sex	Age
7 a	knee	OA	2 / 9	m	75
7 b	knee	OA	2 / 9	m	75
7 d	hip	OA	3 / 9	f	64
7 e	hip	OA	3 / 9	f	64

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.