

pRI 910 DNA (製品コード 3261)

Vector Elements	Location	由来する生物種	実験分類 クラス
ColE1 ori	57...645	<i>Escherichia coli</i>	1
Ri ori	889...5524	<i>Agrobacterium rhizogenes</i>	1
LB (T-DNA left border)	5531...5555	<i>Agrobacterium tumefaciens</i>	1
NOS terminator	5893...6145	<i>Agrobacterium tumefaciens</i>	1
NPTII (neomycin phosphotransferase II)	6419...7213	<i>Escherichia coli</i>	1
NOS promoter	7234...7449	<i>Agrobacterium tumefaciens</i>	1
MCS (multi cloning site)	7620...7679	人工合成核酸	—
RB (T-DNA right border)	7899...7923	<i>Agrobacterium tumefaciens</i>	1
NPTIII (neomycin phosphotransferase III)	8007...8801	<i>Streptococcus faecalis</i>	1