

pRI 201-AN-GUS DNA (製品コード 3266)

Vector Elements	Location	由来する生物種	実験分類 クラス
ColE1 ori	57...645	<i>Escherichia coli</i>	1
Ri ori	889...5524	<i>Agrobacterium rhizogenes</i>	1
LB (T-DNA left border)	5531...5555	<i>Agrobacterium tumefaciens</i>	1
NOS terminator	5893...6145	<i>Agrobacterium tumefaciens</i>	1
NPTII (neomycin phosphotransferase II)	6419...7213	<i>Escherichia coli</i>	1
NOS promoter	7234...7449	<i>Agrobacterium tumefaciens</i>	1
MCS2 (multi cloning site)	7708...7738	人工合成核酸	—
HSP terminator	7744...7992	<i>Arabidopsis thaliana</i>	1
GUS (β -glucuronidase)	8073...9884	<i>Escherichia coli</i>	1
AtADH 5'-UTR	9888...9945	<i>Arabidopsis thaliana</i>	1
35S promoter	9960...10305	Cauliflower Mosaic Virus	1
RB (T-DNA right border)	10951...10975	<i>Agrobacterium tumefaciens</i>	1
NPTIII (neomycin phosphotransferase III)	11059...11853	<i>Streptococcus faecalis</i>	1