

pRI 201-ON-GUS DNA(製品コード 3267)

Vector Elements	Location	由来する生物種	実験分類 クラス
ColE1 ori	57...645	<i>Escherichia coli</i>	1
Ri ori	889...5524	<i>Agrobacterium rhizogenes</i>	1
LB (T-DNA left border)	5531...5555	<i>Agrobacterium tumefaciens</i>	1
NOS terminator	5893...6145	<i>Agrobacterium tumefaciens</i>	1
NPTII (neomycin phosphotransferase II)	6419...7213	<i>Escherichia coli</i>	1
NOS promoter	7234...7449	<i>Agrobacterium tumefaciens</i>	1
MCS2 (multi cloning site)	7708...7738	人工合成核酸	—
HSP terminator	7744...7992	<i>Arabidopsis thaliana</i>	1
GUS (β -glucuronidase)	8073...9884	<i>Escherichia coli</i>	1
OsADH 5'-UTR	9888...9988	<i>Oryza sativa</i>	1
35S promoter	10003...10348	Cauliflower Mosaic Virus	1
RB (T-DNA right border)	10994...11018	<i>Agrobacterium tumefaciens</i>	1
NPTIII (neomycin phosphotransferase III)	11102...11896	<i>Streptococcus faecalis</i>	1