

## Human Cell-Free Protein Expression System (製品コード 3281)

### pT7-IRES Vector

Vector Elements	Location	由来する生物種	実験分類 クラス
Ampr (Ampicillin resistance gene)	201...1061	<i>Escherichia coli</i>	1
ColE1 ori	1232...1820	<i>Escherichia coli</i>	1
PT7	2241...2259	bacteriophage T7	1
EMCV IRES	2271...2822	encephalomyocarditis virus	2
MCS (multi cloning site)	2832...2890	人工合成核酸	—
T7 terminator	2891...2938	bacteriophage T7	1