

**pCold GST DNA(製品コード 3372)**

Vector Elements	Location	由来する生物種	実験分類 クラス
cspA promoter	15..81	<i>Escherichia coli</i>	1
lac operator	82..121	<i>Escherichia coli</i>	1
cspA 5'UTR (untranslated region)	122..255	<i>Escherichia coli</i>	1
TEE (translational enhancing element)	256..270	<i>Escherichia coli</i>	1
His Tag	271..288	人工合成核酸	—
Factor Xa recognition site	289..300	人工合成核酸	—
GST (Glutathione S-transferase) Tag	304..957	<i>Schistosoma japonicum</i>	2
HRV 3C Protease site	964..987	人工合成核酸	—
MCS (multi cloning site)	991..1050	人工合成核酸	—
cspA 3'UTR (untranslated region)	1058..1202	<i>Escherichia coli</i>	1
M13 IG (intergenic region)	1401..1874	M13 phage	1
Ampr (Ampicillin resistance gene)	2061..2921	<i>Escherichia coli</i>	1
ColE1 ori	3081..3695	<i>Escherichia coli</i>	1
lacI	3852..4934	<i>Escherichia coli</i>	1