

B. subtilis Secretory Protein Expression System (製品コード 3380)

pBE-S DNA

Vector Elements	Location	由来する生物種	実験分類 クラス
aprE promoter	7...172	<i>Bacillus subtilis</i>	1
aprE signal peptide	198...284	<i>Bacillus subtilis</i>	1
MCS (multi cloning site)	297...356	人工合成核酸	—
His-Tag	357...374	人工合成核酸	—
ColE1 ori	1698...2297	<i>Escherichia coli</i>	1
Ampr (Ampicillin resistance gene)	2459...3319	<i>Escherichia coli</i>	1
Kanr (Kanamycin resistance gene)	3594...4355	<i>Staphylococcus aureus subsp.</i>	2
pUB ori	4524...5691	<i>Staphylococcus aureus subsp.</i>	2