

**pHEK293 Ultra Expression Vector II (製品コード 3392)**

Vector Elements	Location	由来する生物種	実験分類 クラス
CMV IE promoter	60...567	Human cytomegalovirus	2
TAR (Trans Activation Responsive region)	587...799	HIV-1	3
MCS (multi cloning site)	800...847	人工合成核酸	—
HSV TK polyA	848...947	herpes simplex virus	2
Ampr (Ampicillin resistance gene)	1634...2494	<i>Escherichia coli</i>	1
pUC ori	2654...3268	<i>Escherichia coli</i>	1