

pAUR123 DNA (製品コード 3602)

Vector Elements	Location	由来する生物種	実験分類 クラス
Ampr (Ampicillin resistance gene)	257...1117	<i>Escherichia coli</i>	1
Ori (pBR322 origin)	1273...2063	<i>Escherichia coli</i>	1
CEN4	2893...3003	<i>Saccharomyces cerevisiae</i>	1
ARS1	3382...3437	<i>Saccharomyces cerevisiae</i>	1
AUR1-C	4017...5222	<i>Saccharomyces cerevisiae</i>	1
PADH1	5979...6590	<i>Saccharomyces cerevisiae</i>	1
MCS (multi cloning site)	6591...6624	人工合成核酸	—
TADH1	6621...6816	<i>Saccharomyces cerevisiae</i>	1