

pAUR224 DNA (製品コード 3603)

Vector Elements	Location	由来する生物種	実験分類 クラス
MCS (multi cloning site)	1...75	人工合成核酸	—
SV40 polyA	115...298	Simian Virus 40	2
ars1	301...1506	Schizosaccharomyces pombe	1
ColE1 ori	1875...2529	<i>Escherichia coli</i>	1
Ampr (Ampicillin resistance gene)	2677...3537	<i>Escherichia coli</i>	1
IG	3924...4397	M13 phage	1
aur1R	5003...6271	<i>Schizosaccharomyces pombe</i>	1
PCMV	7101...7608	Human cytomegalovirus	2