

## Tissue Microarray – Normal adult brain tissue

Cat.-No.: 401 1210 Lot: 028P

	a	b	c	d
1	•	•	•	•
2	•	•	•	•
3	•			

Technical Information: 9 spots

- Spot diameter: 1.5 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Code	Tissue of	Sex	Age
1a	H1	Hippocampus/temporal Cortex	female	88
1b	H2	Medulla oblongata	female	88
1c	H3	Cerebrum, white matter	male	72
1d	H4	Cerebellum, white matter	female	88
2a	H5	Pons	female	88
2b	H6	Cerebellum, white matter	male	72
2c	H7	Cerebellum, cortex	male	72
2d	H8	Basal ganglia/internal capsule	female	88
3a	H9	Pons	female	88

H1, H2, H4, H5, H8, H9

Normal tissue from one patient

H3, H6, H7

Normal tissue from a second patient

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.