

## Tissue Microarray – Normal adult cartilage tissue I

Cat.-No.: 401 1221 Lot: 046P

	a	b	c	d	e
1	•	•			
2	•	•	•		
3	•	•	•	•	•
4	•	•	•	•	•
5	•	•	•	•	•

Technical Information: 20 spots

- Spot diameter: 2.0 mm
- Fixation in 4 % paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Localisation	Sex	Age
1a	hip	female	50
1b	hip	female	50
2a	hip	male	21
2b	hip	male	21
2c	trachea	male	56
3a	trachea	male	56
3b	knee	male	15
3c	knee	male	15
3d	trachea	male	25
3e	trachea	male	25
4a	trachea	female	53
4b	trachea	female	53
4c	hip	female	53
4d	hip	female	53
4e	trachea	female	57
5a	trachea	female	57
5b	trachea	male	72
5c	trachea	male	72
5d	knee	male	76
5e	knee	male	76

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.