

## Tissue Microarray - Human Cervical Carcinoma

Cat.-No.: 401 2203 Lot: 017P

	a	b	c	d	e	f
1	•	•	•	•	•	•
2	•	•	•	•	•	•
3	•	•	•	•	•	•
4	•	•	•	•	•	•
5	•	•	•	•	•	•
6	•	•	•	•	•	•
7	•	•	•	•	•	•
8						
9	•	•	•	•	•	
10	•	•	•	•	•	

Technical Information: 52 spots

- Spot diameter: 2 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnostics	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
1a	SCC	female	48	4	0	X	3
1b	SCC	female	30	2b	0	X	3
1c	SCC	female	31	X	1	1	3
1d	SCC	female	39	1b	0	X	3
1e	ADC	female	46	2b	0	X	2
1f	SCC	female	32	4	0	X	3
2a	SCC	female	69	2	X	X	2
2b	SCC	female	45	1b	0	X	2
2c	SCC	female	61	1b	0	X	2
2d	SCC	female	41	1b	0	X	2
2e	SCC	female	44	1b	0	X	2
2f	SCC	female	52	2b	1	X	3
3a	SCC	female	43	1b	0	X	2
3b	SCC	female	43	2b	0	X	2
3c	SCC	female	38	y1b	0	X	3
3d	SCC	female	36	2b	1	X	2
3e	SCC	female	62	2	0	X	2
3f	SCC	female	25	1b	1	X	3
4a	SCC	female	47	2b	1	1	2
4b	lymph node metastasis	female	47	2b	1	1	2
4c	SCC	female	72	2b	X	X	3
4d	SCC	female	40	1b	0	X	2
4e	SCC	female	38	1 a1	X	x	2
4f	SCC	female	68	1 b1	0	x	2

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnostics	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
5a	SCC	female	50	2b	1	X	3
5b	lymph node metastasis	female	50	2b	1	X	3
5c	SCC	female	43	1b	1	X	2
5d	SCC	female	69	1	X	X	3

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnostics	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
5e	SCC	female	56	2b	0	X	3
5f	SCC	female	32	2b	1	X	2
6a	SCC	female	46	2a	0	X	3
6b	SCC	female	70	2b	0	X	2
6c	SCC	female	62	X	X	1	3
6d	Peritumorous inflammatory tissue	female	49	y2b	X	1	3
6e	SCC	female	35	1b2	0	X	2
6f	SCC	female	69	4	0	X	2
7a	SCC	female	42	2b	1	X	2
7b	SCC	female	58	2b	X	X	3
7c	SCC	female	44	1b1	X	X	2
7d	SCC	female	33	1b1	0	X	2
7e	SCC	female	37	1b1	0	X	3
7f	SCC	female	59	1bi	0	X	3
9a	CIN III lesion	female	50				
9b	CIN III lesion	female	63				
9c	Cervix mucosa near CIN III lesion	female	40				
9d	CIN III lesion	female	29				
9e	Cervix mucosa near CIN III lesion	female	26				
10a	normal mucosa	female	30				
10b	normal mucosa	female	32				
10c	normal mucosa	female	35				
10d	normal cervix, fibrous tissue	female	36				
10e	normal mucosa	female	35				

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.