

Tissue Microarray - Human Autoimmune Tissue

Cat.-No.: 401 3301 Lot: 134P

Slide Label						
	a	b	c	d	e	f
1	•	•	•	•	•	•
2	•	•	•	•	•	•
3	•	•	•	•		
4	•	•	•	•		
5	•	•	•	•		
6	•	•	•			
7	•	•	•	•	•	•
8	•	•	•	•	•	•
9	•	•	•	•		

Technical Information: 43 Spots

- Spot diameter: 2.0 mm
- Fixation in 4% paraformaldehyde in PBS
- Paraffin embedded

Tissue type validated by immunohistochemistry

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Tissue	Sex	Age
1a	Hashimoto-Thyreoiditis	Thyroid	f	37
1b	Hashimoto-Thyreoiditis	Thyroid	f	37
1c	Hashimoto-Thyreoiditis	Thyroid	f	13
1d	Hashimoto-Thyreoiditis	Thyroid	f	13
1e	Hashimoto-Thyreoiditis	Thyroid	f	59
1f	Hashimoto-Thyreoiditis	Thyroid	f	59
2a	Primäre biliäre Zirrhose	Liver	f	47
2b	Primäre biliäre Zirrhose	Liver	f	47
2c	Sjögren	Submandibularis	f	72
2d	Sjögren	Submandibularis	f	72
2e	Sjögren	Labium oris	m	70
2f	Sjögren	Labium oris	m	70
3a	Sinusitis & Eosinophilie	Nasal mucosa	m	78
3b	Sinusitis & Eosinophilie	Nasal mucosa	m	78
3c	Sinusitis & Eosinophilie	Nasal mucosa	f	35
3d	Sinusitis & Eosinophilie	Nasal mucosa	f	35
4a	Rheumatoide Arthritis	Synovial	f	46
4b	Rheumatoide Arthritis	Synovial	f	46
4c	Psoriasis	Synovial	m	21
4d	Psoriasis	Synovial	m	37
5a	Morbus Crohn	Sigma	f	27
5b	Morbus Crohn	Sigma	f	27
5c	Sarkoidose	Lung	m	47
5d	Sarkoidose	Lung	m	47

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Tissue	Sex	Age
6a	Granulomat. Thyreoiditis de Quervain	Thyroid	m	40
6b	Granulomat. Thyreoiditis de Quervain	Thyroid	m	40
6c	Lupus erythematodes	Dermis	m	65
7a	normal	Thyroid	m	53
7b	normal	Thyroid	m	53
7c	normal	Liver	m	61
7d	normal	Liver	m	61
7e	normal	Salivary gland	f	51
7f	normal	Salivary gland	f	51
8a	normal	Dermis	f	36
8b	normal	Dermis	f	36
8c	normal (Sinusitis w/o. Eosinophilie)	Nasal mucosa	m	36
8d	normal (Sinusitis w/o. Eosinophilie)	Nasal mucosa	m	36
8e	normal	Sigma	m	58
8f	normal	Sigma	m	58
9a	normal	Lung	m	47
9b	normal	Lung	m	47
9c	normal	Synovial	f	62
9d	normal	Synovial	f	62

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.