

pAAV-ZsGreen1 Vector (製品コード 6231)

pAAV-ZsGreen1 Vector

Vector Elements	Location	由来する生物種	実験分類 クラス
5'ITR (5' Inverted terminal repeat)	237..377	Adeno Associated Virus	1
CMV Promoter	446..953	Human cytomegalovirus	2
Human b-globin intron	988..1480	Homo sapiens	1
ZsGreen1	1499..2194	Zoanthus sp.	1
Human growth hormone(hGH) polyA signal	2204...2682	Homo sapiens	1
3' ITR (3' Inverted terminal repeat)	2695...2835	Adeno Associated Virus	1
Ampr (Ampicillin resistance gene)	3811...4776	<i>Escherichia coli</i>	1
pUC ori (pUC origin of replication)	4947...5535	<i>Escherichia coli</i>	1