

pAAV1-ZsGreen1 One-Vector (製品コード 6263)

Vector Elements	Location	由来する生物種	実験分類 クラス
AAV2 Rep	131..1996	Adeno Associated Virus	1
AAV1 Cap	2013..4223	Adeno Associated Virus	1
Adenovirus E2A	4324..9659	Adenovirus Type5	2
Adenovirus E4	9659..12860	Adenovirus Type5	2
Adenovirus VA	12869..13599	Adenovirus Type5	2
5'ITR	13685..13825	Adeno Associated Virus	1
CMV promoter	13894..14401	Human cytomegalovirus	2
Human beta-globin intron	14436..14928	Homo sapiens	1
ZsGreen1	14976..15671	Adeno Associated Virus	1
Human growth hormone(hGH) polyA signal	15691..16167	Homo sapiens	1
3'ITR	16180..16320	Adeno Associated Virus	1
KmR(kanamycin resistance gene)	17180..18101	Escherichia coli	1
pUC ori	18272..18860	Escherichia coli	1