

pAAV5-ZsGreen1 One-Vector (製品コード 6265)

Vector Elements	Location	由来する生物種	実験分類 クラス
AAV2 Rep	131..1996	Adeno Associated Virus	1
AAV5 Cap	2013..4187	Adeno Associated Virus	1
Adenovirus E2A	4288..9623	Adenovirus Type5	2
Adenovirus E4	9623..12824	Adenovirus Type5	2
Adenovirus VA	12833..13563	Adenovirus Type5	2
5'ITR	13649..13789	Adeno Associated Virus	1
CMV promoter	13858..14365	Human cytomegalovirus	2
Human beta-globin intron	14400..14892	Homo sapiens	1
ZsGreen1	14940..15635	Adeno Associated Virus	1
Human growth hormone(hGH) polyA signal	15655..16131	Homo sapiens	1
3'ITR	16144..16284	Adeno Associated Virus	1
KmR(kanamycin resistance gene)	17144..18065	Escherichia coli	1
pUC ori	18236..18824	Escherichia coli	1