

pAAV8-ZsGreen1 One-Vector (製品コード 6268)

Vector Elements	Location	由来する生物種	実験分類 クラス
AAV2 Rep	131..1996	Adeno Associated Virus	1
AAV8 Cap	2013..4229	Adeno Associated Virus	1
Adenovirus E2A	4330..9665	Adenovirus Type5	2
Adenovirus E4	9665..12866	Adenovirus Type5	2
Adenovirus VA	12875..13605	Adenovirus Type5	2
5'ITR	13691..13831	Adeno Associated Virus	1
CMV promoter	13900..14407	Human cytomegalovirus	2
Human beta-globin intron	14442..14934	Homo sapiens	1
ZsGreen1	14982..15677	Adeno Associated Virus	1
Human growth hormone(hGH) polyA signal	15697..16173	Homo sapiens	1
3'ITR	16186..16326	Adeno Associated Virus	1
KmR(kanamycin resistance gene)	17186..18107	Escherichia coli	1
pUC ori	18278..18866	Escherichia coli	1