



Repli-tech Co., Ltd.
<https://repli-tech.co.jp/>
techinfo@repli-tech.co.jp

Anti-Sendai virus Rabbit pAb IgG

PRODUCT SPECIFICATION

Product Name	Anti-Sendai virus Rabbit pAb IgG
Product Number	A-0100
Clonality	Polyclonal
Isotype	IgG
Host	Rabbit
Antigen	Sendai virus strain Z
Purification	Antibody purification using Protein A
Reactivity	Sendai virus
Applications	Immunofluorescence (1:500), ELISA
Buffer	50% glycerol and PBS. 0.025% sodium azide is added as preservative.
Size	100 μ L
Storage	Store at -20°C or - 80°C. Avoid freeze / thaw cycle. Valid for 12 months from date of delivery.

For research use only. Not intended for pharmaceutical development, diagnostic, therapeutic or any in vivo use. No products from Repli-tech may be resold, modified for resale or used to manufacture commercial products without prior written approval from Repli-tech Co., Ltd.