

## Our commitment to quality

### ISO 18385 Forensic DNA Grade Certificate

The ISO 18385:2016 “Minimizing the risk of human DNA contamination in products used to collect, store and analyze biological material for forensic purposes – requirements” applies to the manufacturing of consumables for forensic DNA purposes including DNA extraction kits. The ultimate aim is to minimize the risk of human DNA contamination in forensic workflows.

MACHEREY-NAGEL has committed to these standards for quality by implementing all requirements of this international standard. In addition we also meet, where feasible, the guidelines of the European Network of Forensic Science Institutes (ENFSI), the Scientific Working Group on DNA Analysis Methods (SWGDM\*) and the Biology Specialist Advisory Group (BSAG). Key elements of the manufacturing comprise:

- **Production of forensic products:**  
Production processes designed to minimize the risks of human DNA contaminations. High degree of automation and dedicated, strictly controlled work environments including positive air pressure for critical manufacturing processes. Strict adherence to cleaning procedures, SOPs and extensive training of staff.
- **Staff elimination database:**  
For exclusion of contamination, MN maintains a staff elimination database containing DNA profiles of all employees involved in the production of forensic products.
- **Forensic quality management system:**  
For highest forensic quality, standards are continuously updated by a dedicated MN Forensic team. Corrective and preventive actions guarantee the highest forensic standards.
- **Post production:**  
All consumables used in our forensic product line are treated (post-production) with ethylene oxide (EO), to guarantee the absence of detectable human DNA. Validated batch-release testing includes tests for the absence of human DNA according to ISO 18385 Annex A.

By product labeling with the “ISO 18385 Certified forensic quality” logo, MACHEREY-NAGEL declares the ISO 18385:2016 compliant manufacturing of these products.

\*SWGDM Contamination Prevention and Detection Guidelines for Forensic DNA Laboratories, 01/12/2017

Düren, July 11, 2022

MACHEREY-NAGEL GmbH & Co. KG



Management System  
EN ISO 13485:2016  
ISO 9001:2015  
  
www.tuv.com  
ID 9000008338

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany · Tel.: +49 24 21 969-0 · E-Mail: info@mn-net.com · www.mn-net.com

Amtsgericht Düren, HRA 1164

GmbH eingetragen beim Amtsgericht Düren, HRB 526

Persönlich haftende Gesellschafterin:

MACHEREY-NAGEL GmbH, Düren

Geschäftsführerin: C. Wagner

Ust.-IdNr.: DE 122277751

Banken: Deutsche Bank AG, Düren

(BLZ 395 700 61) 8020786

Santander Bank, Frankfurt

(BLZ 500 333 00) 9100400206

Commerzbank Düren

(BLZ 395 400 52) 4370300

Sparkasse Düren

(BLZ 395 501 10) 42300608

(IBAN DE66 3957 0061 0802 0786 00) (BIC DEUTDE3395)

(IBAN DE66 5003 3300 9100 4002 06) (BIC SCFBDE33XXX)

(IBAN DE65 3954 0052 0437 0300 00) (BIC COBADEFFXXX)

(IBAN DE56 3955 0110 0042 3006 08) (BIC SDUEDE33XXX)